

SUPERCLOSER SC410

POLYBAG-IN-BOX DECUFFER FOLDER & TOP TAPER



The Supercloser SC410 mechanically un-cuffs a bag in a box clear of the case flaps, opens the side flaps, folds the bag into a neat two-way fold to protect the filled product, closes the flaps and tapes the top of the cases.





TECHNICAL SPECIFICATIONS SUPERCLOSER SC410:

Case range:	250mm x 200mm x 150mm (L x W x H) min. 500mm x 400mm x 500mm (L x W x H) max.
Speed:	Up to 12 cases per min.
Tape Applicator:	OK-55 high quality taping head with open tape feed path, suitable for all standard 50mm wide pressure sensitive tapes up to 350mm roll diameter.
Bag materials:	All standard polyethylene films (HDPE, LDPE, LLDPE, co-extrusions)
Bag gauges:	25-100 micron depending on material.
Pneumatic:	6-8 bar, 1/2" supply.
Finish:	Stainless steel frame, covers, contact parts and electrical enclosure. Internal components stainless steel, aluminum or plated.
Machine handling:	Right to left or left to right to requirement.

SUPERCLOSER SC410 OPTIONS:

- OK-75 taping head for 75mm wide tape.
- Bag top heat seal assembly.
- Top sealing with hot-melt glue.
- Smaller or larger case sizes.

FEATURES

Stainless steel finish.	Clean and durable, maintains hygiene.
Heavy duty conveyor system.	Highest speed at low maintenance.
In-feed gate.	Separates cases for automatic operation.
Opens the side flaps.	More positive folding process.
Fully interlocking guarding.	Operator safety and OSHA compliance.
Touch screen display.	Operator friendly and easy troubleshooting

BENEFITS



Our machines can be used to pack almost any product:

Fruits, nuts, vegetables, scallions, meats, poultry, seafood, cereals, crackers, rice, pasta, tortillas, snack food, sugar, chocolate, herbs, spices, seed, coffee, tea, bakery goods, gypsum, asphalt, caps, closures, plastics, screws, medical devices, cutlery, butter, powders, pharmaceutical devices, biomedical devices, pellets, powders, crack filler, margarine, shortening, oils and fats, dough, candy.

